



Acquisition guide



Transmural Biotech S.L. ; Calle Juan Ignacio Luca de Tena 12, 28027 Madrid España




DATE: 28/01/2026

NUMBER OF EDITION: 2


Acquisition guide for thorax RX



- Visualization of the full lung fields.
- The image should fully include both apex and costodiaphragmatic sinuses.
- Preferably Posterior-Anterior (PA) projection than Anterior-Posterior (AP) to avoid heart shadow.



Lateral projection not allowed.



Pacemakers or other implants in the thoracic area not allowed.